

# CARLO Aiu, Classic 185




**beka** DIFFERENT  
BY DESIGN





Watch product  
video!



A safe and reliable  
passive lift for your residents.

Our mobile sling lift CARLO Alu, Classic 185, provides your care staff with optimum support for all daily lifting and transfer routines and enables them to work ergonomically, safely and efficiently at all times without placing additional strain on the musculoskeletal system by manually lifting and repositioning the resident.

The extending arm with 4-point loop bar can be lowered to the floor using the hand control, ensuring that a person who has fallen from the floor can be picked up safely and easily by just one carer. The chassis can be spread electrically and therefore offers a high level of operating comfort. Due to the low height of the chassis, the CARLO Alu, Classic 185 model, can easily be used to drive under both bathtubs and care beds.



Group 3

Group 4

Group 5





## Stable 4-point loop bar

The CARLO Alu, Classic 185 belt bar is suitable for all loop belts and has a safety locking mechanism for the belt holder. The distance between the mounting bar and the lifting pole is large so that residents have more space during the transfer.



## Anti-rotation mechanism

The anti-rotation mechanism (ARM) prevents unpleasant swinging or rotating movements during lifting and transfer.



## Ideal for bathing and showering

The passive lift can be used in the bathroom without restriction by using high-quality aluminum, both for showering and bathing your residents. The extremely lightweight design makes it easy to maneuver without any effort. You can purchase the appropriate (bath) straps from us.

### Technichal Data

Chassis ( W x H )	1.400 mm x 125 mm
Width inside ( closed/spread )	540 mm / 1.040 mm
Width outside ( closed/spread )	700 mm / 1.180 mm
Stroke	1.400 mm
SWL	185 kg

### Product information

Stable aluminum-steel construction
Hand control with battery charge status display
4-point loop handle with safety locking mechanism
Electric height adjustment and electrically spreadable chassis
Anti-rotation mechanism (ARM)

### Available Options



Extra battery 24 V



Wall charging unit



Digital patient scale



Low chassis

### Clip Slings



Transfer slings  
With headrest  
Padded leg loops  
High quality nylon material



Bathing slings  
With headrest  
Padded leg loops  
High quality nylon-mesh material



Toilet slings  
With headrest  
Padded leg loops  
Arm padding



Soft slings  
With headrest  
Padded leg loops  
Breathable nylon/cotton blend



**BEKA Hospitec GmbH** | Am Rübenmorgen 3 | 35582 Wetzlar-Dutenhofen  
Tel +49 641 922 20-0 | E-Mail [info@beka-hospitec.com](mailto:info@beka-hospitec.com)  
[www.beka-hospitec.com](http://www.beka-hospitec.com)

